



Sulphuric acid standard solution, 0.020 N, 1 L

Product #:

20353

Hazardous

AED Price:
Available

Contact Hach



Sulfuric Acid Standard Solution, 0.020 N (N/50), 1 L

Specifications

Application:	Remove Chlorine Interference
Concentration:	0.020 N
Description:	Remove Interfering substance
EPA compliant:	Yes
Formula:	H ₂ SO ₄ 0.02 N
Measuring range:	0 - 10
Measuring range (2):	0 - 5000 mg/L CaCO ₃
Method number:	8288 (Acid)
Method number :	8289 (Base)
Package Type:	Bottle
Parameter:	Alkalinity
Platform:	Solution
Quantity:	1 L

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)