



## QAC reagent 1, 0.2- 5.0 mg/L CTAB, 50 powder pillows

<b>Product #:</b>	<b>2401066</b>
AED Price:	Contact Hach
Available	

Powder reagent in sealed package for the determination of Quaternary Ammonium Compounds (QAC).

---

### Specifications

Method:	Direct Binary Complex
Method number:	8336
Number of tests:	50
Parameter:	Quaternary Ammonium Compounds
Platform:	Powder Pillows
Quantity:	50/pk

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)