



QAC reagent 2, 0.2- 5.0 mg/L CTAB, 25 powder pillows

Product #:

2401268

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for the determination of Quaternary Ammonium Compounds (QAC).

Specifications

Method:	Direct Binary Complex
Method number:	8337
Number of tests:	12.5
Parameter:	Quaternary Ammonium Compounds
Platform:	Powder Pillows
Quantity:	25/pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)