



Hexamethylenetetramine buffer reagent, 100 powder pillows

Product #:

2603999

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for the determination of Aluminium with ECR Method. For 20 mL samples.

Specifications

Measuring range:	0.002 - 0.250 Al
Method:	Eriochrome Cyanine R
Method number:	8326
Number of tests:	100
Parameter:	Aluminum
Platform:	Powder Pillows
Quantity:	100 pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)