



Potassium 2 Reagent Solution Pillows, pk/25

Product #:

1432298

AED Price:
Available

Contact Hach

Hazardous



Solution pillows in sealed packages for the analysis of Potassium, 0.1-7.0 mg/L K. For 25mL sample.

Specifications

Method:	Tetraphenylborate
Method number:	8049
Number of tests:	25
Parameter:	Parameter
Quantity:	25 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)