



## Potassium 3 Reagent Powder Pillows, pk/100

**Product #:**

**1432399**

**Hazardous**

AED Price:  
Available

Contact Hach



Powder reagent in sealed package for the analysis of Potassium, 0.1-7.0mg/l K.

For potassium determination in water, soil and plant tissues by the tetraphenylborate method. Hach Method 8049 Sample size: 25 mL. Pack of 100 powder pillows.

---

### Specifications

Method:	Tetraphenylborate
Number of tests:	100
Platform:	Powder Pillows
Quantity:	100 /pk

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)