



Potassium Reagent Set, 0.1-7.0mg/L K

Product #:

2459100

Hazardous

AED Price:
Available

Contact Hach



Powder reagents and reagent solution in sealed packages for determination of Potassium by the Tetraphenylborate method. For 25mL samples; pk/100.

Specifications

Description:	Potassium Reagent Set
Instrument:	All except DR/800, PC II
Measuring range:	0.1 - 7.0 mg/L K
Method:	Tetraphenylborate
Method number:	8049
Number of tests:	100
Parameter:	Potassium
Platform:	Powder Pillows
Special remarks:	Set contains: Potassium 1 Reagent Powder Pillows (1432198), Potassium 2 Reagent Solution Pillows (1432298), Potassium 3 Reagent Powder Pillows (1432399)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)