



Sodium Hydroxide Digital Titrator Cartridge, 0.3636 N

Product #:

1437801

Hazardous

AED Price:
Available

Contact Hach



Ready-to-use titrant solution for the determination of Carbon Dioxide, 10-100mg/l CO₂.

For carbon dioxide determination by Digital Titration

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Measuring range:	10 - 100 mg/L CO ₂
Method:	Phenolphthalein
Method number:	8205
Number of tests:	40 - 100
Parameter:	Carbon Dioxide
Platform:	Digital Titrator
Titrant:	0.3636 N NaOH

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)