



Sodium Hydroxide Digital Titrator Cartridge, 1.600 N

Product #:

1437901

Hazardous

AED Price:
Available

Contact Hach



Ready-to-use titrant solution for Acid-Base Titration (1-10 meq/l acid) and determination of Acidity (100-4000mg/l CaCO_3).

For acidity determination by Digital Titration, 100 - 4000 mg/L CaCO_3

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

EPA compliant:	Yes
Measuring range:	1 - 4000 mg/L CaCO_3
Measuring range (2):	1 - 10 meq/L
Method:	Methyl Orange
Method number:	8200 (ácido)
Method number :	8202 (P)
Number of tests:	40 - 100
Parameter:	Acidity
Platform:	Digital Titrator
Titrant:	1.600 N NaOH

Required Accessories

- Protective gloves size L, blue, nitrile, powder-free, 50 pieces (Item SM743L)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)