



Sulfuric acid cartridge for digital titrator, 0.1600 N

Product #:

1438801

Hazardous

AED Price:
Available

Contact Hach



Ready-to-use titrant solution for the determination of alkalinity, 10 - 160 mg/L CaCO_3

Specifications

EPA compliant:	Yes
Measuring range:	10 - 160 mg/L CaCO_3
Method:	Phenolphthalein
Method number:	8203
Number of tests:	40 - 100
Parameter:	Alkalinity, Total
Platform:	Digital Titrator
Titrant:	0.1600 N H_2SO_4

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)