



Hexamethylenetetramine Reagent Powder Pillows, pk/50

Product #:

187866

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for the determination of Aluminium, ECR Method, 0-0.220 mg/L Al.

Specifications

Instrument:	DR2700/DR2800/DR3800/DR3900/DR5000/DR6000
Method:	Eriochrome Cyanine R
Method number:	8326
Parameter:	Aluminum
Platform:	Powder Pillows
Quantity:	50 pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)