



Porphyrim 1 Reagent Powder Pillows

Product #:	2603549
AED Price:	Contact Hach
Available	

Powder reagent in sealed package for the determination of Copper with Porphyrim method, 2-210µg/L Cu.

Specifications

Method:	Porphyrim
Method number:	8143
Number of tests:	50 - 100
Parameter:	Copper
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Copper Masking Reagent Powder Pillows, 10 mL, pk/100 (Item 2603449)
- Porphyrim 2 Reagent Powder Pillows (Item 2603649)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)