



Chlorine (Free and Total) Reagent Set, 0-3.5 mg/L Cl_2

Product #:

2438800

Hazardous

AED Price:
Available

Contact Hach



For determination of Free and Total Chlorine by the DPD method, using visual test kits. Set contains: DPD Free Chlorine PP, pk/100 and DPD Total Chlorine PP pk/100. Range: 0 - 3.5 mg/L Cl_2 . Sample Size: 5 mL.

Boric Acid free!

Specifications

Measuring range:	0 - 3.5 mg/L Cl_2
Method:	DPD
Number of tests:	100
Sample size:	5 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)