



Chlorine dioxide reagent set 0.02 - 0.50 mg/L ClO_2 , 100 pcs

Product #: HPT240
AED Price: Contact Hach

Amaranth Method

Specifications

Description:	liquid reagent set for Chlorine dioxide
Measuring range:	0.02 - 0.50 mg/L ClO_2
Method:	Amaranth Method
Number of tests:	100
Parameter:	Chlorine dioxide
Special remarks:	Calibrated for DR 820/50/90, 2000/10, 4000, 2500 and POCKET Colorimeter Chlorine
Storage conditions:	15 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)