



Chlorine, liquid DPD reagent set, 0.02-2.00mg/L Cl_2

Product #:

HPT310

Hazardous

AED Price:

Contact Hach



Specifications

Description:	Free / total Chlorine
Measuring range:	0.02 - 2.00 mg/L Cl_2
Method:	DPD
Method number:	RS
Number of tests:	100
Parameter:	Chlorine
Platform:	Solution
Special remarks:	for the evaluation 1 inch round vials are needed
Standard method:	ISO 7393-2
Storage conditions:	15 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)