



## Swiftest DPD Free chlorine reagent dispenser and reagent vial

**Product #:**

**2802300**

**Hazardous**

AED Price:  
Available

Contact Hach



For the determination of Free Chlorine by the DPD method, 0.02-2.00 mg/L Cl<sub>2</sub>.

---

### Specifications

Instrument:	All except DR/3000, DR/2000
Measuring range:	0.02 - 2.00 mg/L Cl <sub>2</sub>
Method:	DPD
Method number:	8021
Number of tests:	Approximately 250
Parameter:	Chlorine, free
Platform:	Swiftest

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)