



Chlorine Dioxide Reagent Set, 0.01-1.00 mg/L ClO_2

Product #:

2242300

Hazardous

AED Price:
Available

Contact Hach



For the determination of Chlorine Dioxide, Chlorophenol Red Method. For 50mL sample. 100 tests.

Specifications

Instrument:	All except DR/800 Series, DR 900, PC II
Measuring range:	0.01 - 1.00 mg/L ClO_2
Method:	Chlorophenol Red
Method number:	8065
Number of tests:	100
Parameter:	Chlorine dioxide

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)