



Free Chlorine Reagent, Accuvac Ampules, DPD, pk/25

Product #:

2502025

Hazardous

AED Price:
Available

Contact Hach



Ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Free Chlorine, 0.02-2.00 mg/L Cl₂.

Specifications

EPA compliant:	Yes
Measuring range:	0.02 - 2.00 mg/L Cl ₂
Method:	DPD
Method number:	8021
Method number :	10126
Number of tests:	25
Parameter:	Chlorine, free, Chlorine dioxide
Platform:	Accuvac
Quantity:	25 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)