



Alkali solution for colorimetric hardness, 100 mL

Product #:

2241732

Hazardous

AED Price:
Available

Contact Hach



For colorimetric Calcium and Magnesium determination, 0.05-4.00 mg/L Ca as CaCO₃, 0.05-4.00 mg/L Mg as CaCO₃.

Specifications

Method:	Calmagite Colorimetric
Method number:	8030
Number of tests:	100
Parameter:	Hardness, Ca and Mg
Platform:	Solution
Quantity:	100 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)