



Calcium / Magnesium indicator solution, 100 mL

Product #:

2241832

Hazardous

AED Price:
Available

Contact Hach



Used for colorimetric Calcium and Magnesium determinations by the Calmagite colorimetric method.

Specifications

Method:	Calmagite Colorimetric
Method number:	8030
Number of tests:	100
Parameter:	Calcium and Magnesium Hardness
Platform:	Solution
Quantity:	100 mL

Required Accessories

- Alkali solution for colorimetric hardness, 100 mL (Item 2241732)
- EDTA solution, 1M, 50mL SCDB (Item 2241926)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)