



Hardness Reagent Set, Calmagite Colorimetric Method

Product #:

2319900

Hazardous

AED Price:
Available

Contact Hach



For the determination of Hardness, 0.05-4.00mg/L Ca as CaCO₃/0.05-4.00mg/L Mg AS CaCO₃. 100 tests

Specifications

Instrument:	All except PC II
Measuring range:	0.05 - 4.00 mg/L CaCO ₃
Measuring range (2):	0.05 - 4.00 mg/L Mg as MgCO ₃
Method:	Calmagite Colorimetric
Method number:	8030
Number of tests:	100
Parameter:	Hardness, Ca and Mg
Platform:	Solution
Special Features:	Reagent set

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)