



Silica Reagent Set, 1-100 mg/L SiO₂

Product #:

2244300

Hazardous

AED Price:
Available

Contact Hach



Powder reagents in sealed packages for the determination of Silica, Silicomolybdate Method. For 25 mL sample. 100 tests.

Specifications

Instrument:	DR/3000, DR/2000, DR/800
Measuring range:	1 - 100 mg/L SiO ₂
Method:	Silicomolybdate
Method number:	8185
Number of tests:	100
Parameter:	Silica
Platform:	Powder Pillows
Sample size:	25 mL
Special Features:	HR, 25 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)