



Silica Reagent Set, Bulk Solution, 3-1000 µg/L SiO₂

Product #:

2553500

Hazardous

AED Price:
Available

Contact Hach



For the determination of ultra low range Silica, Heteropoly Blue Method. Approximately 100 tests.

Specifications

Instrument:	All except DR 900, DR/800 Series, PC II
Measuring range:	3 - 1000 µg/L SiO ₂
Method:	Heteropoly Blue
Method number:	8282
Number of tests:	100
Parameter:	Silica
Platform:	Solution
Sample size:	Bulk
Special Features:	UL, Bulk, 100 tests

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)