



Silica Reagent Set, 3-1000 µg/L SiO₂

Product #:

2581400

Hazardous

AED Price:

Contact Hach

Call for ship date



For ultralow range Silica determination by the Heteropoly Blue method.

Specifications

Instrument:	All except DR/3000, DR/800, PC II
Measuring range:	3 - 1000 µg/L SiO ₂
Method:	Heteropoly Blue
Method number:	8282
Number of tests:	40
Parameter:	Silica
Platform:	Solution
Sample size:	25 mL
Special Features:	UL, 40 tests

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)