



HydraVer2 Hydrazine reagent solution, 4 - 600 µg/L, 100 mL

Product #:

179032

Hazardous

AED Price:
Available

Contact Hach



For the determination of Hydrazine, p-Dimethylaminobenzaldehyde method. For 10 mL sample. 100 mL bottle.

Specifications

Description:	HydraVer 2 Hydrazine Reagent
Instrument:	All except DR/850, DR/820, PC II
Measuring range:	4 - 600 µg/L N ₂ H ₄
Method:	p-Dimethylaminobenzaldehyde
Method number:	8141
Number of tests:	100
Parameter:	Hydrazine
Platform:	Solution
Quantity:	100 mL
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)