



## COD TNT, Mercury free, 0 - 1500 mg/L O<sub>2</sub>

**Product #:** 2345952  
**AED Price:** Contact Hach

Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

---

### Specifications

Description:	Chemical Oxygen Demand
Instrument:	DR 900
Measuring range:	0 - 1500 mg/L O <sub>2</sub>
Method:	Dichromate without mercury (not approved for USEPA reporting purposes)
Method number:	8000
Number of tests:	25 (1 cuvette for blank needed)
Parameter:	COD
Platform:	Test 'N Tube
Special remarks:	16mm HACH cuvette, not readable at DR2700/ DR2800 / DR3800 / DR5000
Storage conditions:	15 °C - 25 °C (protect from light)

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)