



Dissolved Oxygen Reagent Accuvac Ampules, 6-800 µg/L O₂

Product #:

2501025

Hazardous

AED Price:
Available

Contact Hach



Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Dissolved Oxygen, Indigo Carmine Method.

Specifications

Instrument:	All except DR/820, PC II
Measuring range:	6 - 800 µg/L O ₂
Method:	Indigo Carmine
Method number:	8316
Number of tests:	25
Parameter:	Oxygen, dissolved
Platform:	Accuvac
Quantity:	25 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)