



TOC (Purging Method) TNT, 0.3-20.0 mg/L C

Product #:

2760345

Hazardous

AED Price:
Available

Contact Hach



Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

Specifications

Digestion required:	Yes
Instrument:	DR 900
Measuring range:	0.3 - 20.0 mg/L C
Method:	Direct
Method number:	10129
Number of tests:	25 - 50
Parameter:	TOC
Platform:	Test 'N Tube

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)