



TOC (Purging Method) TNT, 100-700 mg/L C

Product #:

2760445

Hazardous

AED Price:

Contact Hach

Call for ship date



Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

Specifications

Digestion required:	Yes
Instrument:	DR 900
Measuring range:	100 - 700 mg/L C
Method:	Direct
Method number:	10128
Number of tests:	25 - 50
Parameter:	TOC
Platform:	Test 'N Tube

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)