



Total Trihalomethanes Plus Reagent Set (THM)

Product #:

2790800

Hazardous

AED Price:

Contact Hach

No availability shown



Screening tool to optimize treatment processes, to locate areas of high THM levels in distribution systems and to determine levels of trihalomethanes in drinking water samples.

Specifications

Digestion required:	Yes
Instrument:	DR 6000, DR/4000, DR 3900, DR 2800, DR 2700, DR/2400
Measuring range:	10 - 600 µg/L CHCl ₃
Method:	THM Plus
Method number:	10132
Number of tests:	50 - 99
Parameter:	Trihalomethanes

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)