



Mercuric thiocyanate solution (200 mL)

Product #:

2212129

AED Price:
Available

Contact Hach

Hazardous



Reagents

For colorimetric chloride determination by the mercuric thiocyanate method. Range: 0.1 - 25.0 mg/L Cl.

Specifications

Method:	Mercuric Thiocyanate
Method number:	8113
Number of tests:	100
Parameter:	Chloride
Platform:	Solution
Quantity:	200 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)