



Mercuric nitrate cartridge for digital titrator, 2.256N

Product #:

92101

AED Price:
Available

Contact Hach

Hazardous



For titration of chloride, 100 - 8000 mg/L Cl⁻

Specifications

Measuring range:	100 - 8000 mg/L Cl ⁻
Method:	Mercuric Nitrate
Method number:	8206
Number of tests:	40 - 100
Parameter:	Chloride
Platform:	Digital Titrator
Titant:	2.256 N Hg(NO ₃) ₂

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)