



Chloride Reagent Set, 10-8000mg/L Cl

Product #:

2272600

AED Price:
Available

Contact Hach

Hazardous



For determination of chloride by mercuric nitrate Digital Titrator Method (Hach Method 8206). Reagent Set includes Diphenylcarbazone Reagent Powder Pillows, Mercuric Nitrate Digital Titrator Cartridges.

Specifications

Measuring range:	10 - 8000 mg/L Cl ⁻
Method:	Mercuric Nitrate Titration
Method number:	8206
Number of tests:	100
Platform:	Digital Titrator

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)