



Chloride Reagent Set, 10-10000mg/L Cl

Product #:

2288000

AED Price:
Available

Contact Hach

Hazardous



For the determination of chloride by Silver Nitrate Digital Titrator Method (Hach Method 8207). Set contains Chloride 2 Indicator Powder Pillows and two Silver Nitrate Titration Cartridges.

Specifications

Measuring range:	10 - 10000 mg/L Cl ⁻
Method:	Silver Nitrate Titration
Method number:	8207
Number of tests:	100
Platform:	Digital Titrator

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)