



Sodium chloride standard solution, 1990 $\mu\text{S}/\text{cm}$ (NIST), 1 L

Product #:	210553
AED Price:	Contact Hach
Available	

1000 mg/L NaCl. For conductivity probes.

Also suitable for standardising conductivity apparatus, $\mu\text{S}/\text{cm}$ ($\mu\text{mhos}/\text{cm}$).

Specifications

Concentration:	1990 $\mu\text{S}/\text{cm}$
Parameter:	Conductivity
Platform:	Solution
Quantity:	1 L