



Mixed Parameter Quality Control Standard (500 mL)

Product #:	2833149
AED Price:	Contact Hach
Available	

Quality Control Standard (waste water influent inorganics) for the parameters Ammonia, Nitrate, Phosphate, COD, Sulfate and TOC.

The standard solution contains the following concentrations:

- Ammonia: 15 mg/L $\text{NH}_3\text{-N}$
- Nitrate: 10 mg/L $\text{NO}_3\text{-N}$
- COD: 500 mg/L O_2
- Phosphate: 10 mg/L PO_4
- Sulphate: 400 mg/L SO_4
- TOC: 161 mg/L C

Specifications

Concentration:	Nitrate 10 mg/L $\text{NO}_3\text{-N}$
Parameter:	Mixed
Platform:	Solution