



Mixed Parameter Quality Control Standard (500mL)

Product #:	2833249
AED Price:	Contact Hach
Available	

Quality Control Standard (waste water effluent inorganics) for the parameters Ammonia, Nitrate, Phosphate, COD, Sulphate and TOC.

The standard solution contains the following concentrations:

- Ammonia 2.0 mg/L $\text{NH}_3\text{-N}$ / 2.1 mg/L $\text{NH}_4\text{-N}$
- Nitrate 4.0 mg/L $\text{NO}_3\text{-N}$
- Phosphate 2.0 mg/L PO_4
- COD 25 mg/L O_2
- Sulphate 50 mg/L SO_4
- TOC 8 mg/L C

Specifications

Concentration:	COD 25 mg/L O_2
Description:	Wastewater Effluent Inorganics Quality Control Standard, 500 mL
Platform:	Solution
Storage conditions:	10 - 25 °C