



Stablcal Turbidity Standards Calibration Kit, 0 - 7500 NTU

Product #:

2659505

Hazardous

AED Price:

Contact Hach

Ships within 3 days



Six ready to use, stabilised Formazine Primary standard solutions. Kit includes sealed vials of each: <0.1 NTU, 20 NTU, 200 NTU, 1000 NTU, 4000 NTU, 7500 NTU. Calibration Kit, for 2100AN, 2100AN IS, TL2350 and TL2360 IS Laboratory Turbidimeter.

Specifications

Concentration: < 0.1, 20, 200, 1000, 4000, 7500 NTU

Quantity: 6 (set of)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)