



Stablcal Stabilised Turbidity Standard, 10.0 NTU, 500 mL

Product #:

2659949

Hazardous

AED Price:
Available

Contact Hach



Ready to use, stabilised Formazine Primary standard solution in 500 mL bottle, 10 NTU. For TL23, TU52, 2100Q, 2100N, 2100AN series.

Specifications

Concentration:	10 NTU
Description:	STABLCAL Solution 10 NTU
Parameter:	Turbidity
Platform:	Solution
Quantity:	500 mL
Storage conditions:	5 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)