



Stablcal Stabilised Formazine Turbidity Standard, 100 NTU, 500 mL

Product #:

2660249

Hazardous

AED Price:
Available

Contact Hach



Ready to use, stabilised Formazine Primary standard solution, 100 NTU, 500 mL bottle. For TL23, TU52, 2100Q, 2100N, 2100AN series.

Specifications

Concentration:	100 NTU
Platform:	Solution
Quantity:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)