



## Stablcal Stabilised Formazine Turbidity Standard, 1000 NTU, 500 mL

**Product #:**

**2660649**

**Hazardous**

AED Price:  
Available

Contact Hach



Ready to use, stabilised Formazine Primary standard solution in 500 mL bottle, 1000 NTU. For TL23, TU52, 2100Q, 2100N, 2100AN series.

---

### Specifications

|                |          |
|----------------|----------|
| Concentration: | 1000 NTU |
| Platform:      | Solution |
| Quantity:      | 500 mL   |

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)