



## Stablcal Turbidity Standard, 0.50 NTU, 100 mL

**Product #:**

**2698042**

**Hazardous**

AED Price:  
Available

Contact Hach



Ready to use, stabilised Formazine primary standard solution, 0.50 NTU, 100 mL bottle. For TL23, TU52, 2100Q, 2100N, 2100AN series.

### **Precise, pre-mixed Formazin solutions**

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

### **Based on Formazin Primary Standards**

Stablcal Stabilised Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

### **Enhanced stability saves time and labour**

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilised version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

### **Confidence in results**

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

---

## Specifications

Concentration:	0.5 NTU
Platform:	Solution
Quantity:	100 mL

---

## Required Accessories

- Protective gloves size L, blue, nitrile, powder-free, 50 pieces (Item SM743L)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)