



Stablcal Ampule Calibration Kit, Formazin, for 2100Q/QIS portable turbidimeter

Product #:

2971205

AED Price:
Available

Contact Hach

Hazardous



Calibration Standard kit for 2100Q/QIS portable turbidimeter. Ampule kit contains Formazin standards in concentrations: 10, 20, 100, 800 NTU.

Precise, Pre-mixed Formazin Solutions

Based on Formazin Primary Standards

Enhanced Stability Saves Time and Labor

Confidence in results

Specifications

Concentration: 20, 100, 800, 10 NTU (for verification)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)