



Nitrite Reagent Powder Pillows, 2-150 mg/L NO₂

Product #:

221969

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for the determination of Nitrite. 25 mL sample, pk/100.

Specifications

Description:	NitriVer 2 Nitrite Reagent Powder Pillows, 25 mL
Instrument:	DR/2000, DR/3000
Measuring range:	2 - 150 mg/L NO ₂ ⁻
Method:	Ferrous Sulfate
Method number:	8153
Number of tests:	100
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)