



Ammonia Cyanurate Reagent Powder Pillows, 25 mL, pk/50

Product #:

2395566

Hazardous

AED Price:
Available

Contact Hach



Powder reagents in sealed package for the determination of Nitrogen, Ammonia, 0.01-0.50mg/L $\text{NH}_3\text{-N}$. 25mL sample, 50 pcs.

Specifications

Method:	Salicylate
Method number:	8155
Number of tests:	25
Parameter:	Ammonia Nitrogen
Platform:	Powder Pillows
Quantity:	50 /pk