



NitriVer 2 Nitrite Reagent PP, 5 mL, pk/100

Product #:

181369

Hazardous

AED Price:

Contact Hach

Call for ship date



Powder reagent in sealed packages for the determination of high range Nitrite (up to 250 mg/l NO₂).

For determination of high range nitrite by the ferrous sulfate method. Sample Size: 5 mL. pk/100 Powder Pillows.

Specifications

Method:	Ferrous Sulfate
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)