



## Nitrogen, total, TNT, 10-150 mg/L N

**Product #:**

**2714100**

**Hazardous**

**AED Price:**

Contact Hach

**Ships within 3 days**



Cuvette test (Ø 16 mm) completely equipped with all reagents, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube. 25-50 tests

---

### Specifications

Description:	Nitrogen, Total TNT Reagent Set, HR
Digestion required:	Yes
Instrument:	DR 900
Measuring range:	2 - 150 mg/L N
Method:	Persulphate Digestion
Method number:	10072
Number of tests:	25 - 50
Parameter:	Nitrogen, total
Platform:	Test 'N Tube
Special remarks:	Set consists of 2 cartons: 27140-45 and 26721-45  27140-45 = TN Hydroxyde reagent solution (2714045), Persulfate Powder Pillows (2671846)  26721-45 = TN Acid reagent C solution (2672145), Bisulfitreagent A Powder Pillows (2671946), Indicator reagent B Powde
Storage conditions:	10 °C - 25 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)