



## SPADNS 2 Fluoride Reagent, arsenic-free, ACCUVAC, 0.02-2.00 mg/L F

**Product #:**

2527025

Hazardous

AED Price:  
Available

Contact Hach



ACCUVAC ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Fluoride by the SPADNS 2 method. pk/25

---

### Specifications

EPA compliant:	Yes
Measuring range:	0.02 - 2.00 mg/L F <sup>-</sup>
Method:	SPADNS 2
Method Number 1 Chemistry:	10225
Number of tests:	25
Parameter:	Fluoride
Platform:	Accuvac
Quantity:	25 /pk

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)