



Ozone Accuvac, 0.01-1.5 mg/L O₃

Product #:

2518025

Hazardous

AED Price:
Available

Contact Hach



Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of ozone by the Indigo method. pk/25

Specifications

Instrument:	All except DR/820, PC II
Measuring range:	0.01 - 1.50 mg/L O ₃
Method:	Indigo
Method number:	8311
Number of tests:	25
Parameter:	Ozone
Platform:	Accuvac

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)