



## Phosphonate Reagent Set, Powder Pillows, 0.02-125mg/L PO<sub>4</sub>

**Product #:**

**2429700**

AED Price:  
Available

Contact Hach

**Hazardous**



Powder reagents in sealed packages for determination of Phosphonate compounds by UV Persulfate Oxidation Method.

### Specifications

Description:	Phosphonate Reagent Set
Instrument:	All except DR/820
Measuring range:	0.02 - 2.50 mg/L PO <sub>4</sub>
Measuring range (2):	1.0 - 125 mg/L PO <sub>4</sub>
Method:	Persulfate UV Oxidation
Method number:	8007
Number of tests:	100
Parameter:	Phosphonates
Platform:	Powder Pillows
Sample size:	10 mL
Special remarks:	Set contains: Potassium Persulfate Powder Pillows (2084769), PhosVer 3 Powder Pillows (2106069)
Storage conditions:	10 °C - 25 °C

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)